

Protein G Sepharose™ 4 Fast Flow

Protein G Sepharose 4 Fast Flow is an affinity resin with protein G immobilized to Sepharose 4 Fast Flow by the CNBr method.

Protein G binds to the Fc region of IgG from a variety of mammalian species. Protein G Sepharose 4 Fast Flow may be used to isolate and purify classes, subclasses and fragments of immunoglobulins from any biological fluid or cell culture medium. Protein G Sepharose 4 Fast Flow is extremely useful for isolation of immune complexes. See [Section 5 Immunoprecipitation, on page 11](#) for detailed instructions on immunoprecipitation application.

The potential applications of protein G include practically all of the current and projected applications of protein A. Protein G and protein A, however, have different IgG binding specificities, dependent on the origin of the IgG. Compared to protein A, protein G binds more strongly to polyclonal IgG, for example, from cow, sheep and horse. Furthermore, unlike protein A, protein G binds polyclonal rat IgG, human IgG3 and mouse IgG1.



Table of Contents

1	BioProcess™ resins	3
2	Product description	3
3	Operation	6
4	Cleaning-in-place and Sanitization	10
5	Immunoprecipitation	11
6	Storage	16
7	Ordering information	17

Please read these instructions carefully before using the products.

Safety

For use and handling of the products in a safe way, please refer to the Safety Data Sheets.

1 BioProcess™ resins

BioProcess resins are developed and supported for production scale chromatography. All BioProcess resins are produced with validated methods and are tested to meet manufacturing requirements. Secure ordering and delivery routines give a reliable supply of resins for production scale. Regulatory Support Files (RSF) are available to assist process validation and submissions to regulatory authorities. BioProcess resins cover all purification steps from capture to polishing.

2 Product description

Recombinant protein G, M_r 17 000, from GE is produced in *E. coli* and contains two IgG binding regions. The albumin binding region of native protein G has been genetically deleted, thereby avoiding undesirable cross-reactions with albumin. The pI of protein G is 4.1 and the pH stability 2 to 10.

Table 1. Characteristics of Protein G Sepharose 4 Fast Flow

Matrix	Highly cross-linked spherical agarose, 4%
Ligand density	Approx. 2 protein G/mL drained resin
Dynamic binding capacity ¹	Approx. 18 mg human IgG/mL drained resin
Average particle size	90 µm (45 to 165 µm)
Recommended flow rate ²	50 to 300 cm/h
Max. operating pressure	0.1 MPa (1 bar, 14 psi)
Chemical stability	Stable in: 7 days at 37°C: 1 M acetic acid, pH 2.0; 20 mM sodium phosphate; 1% SDS, pH 7.0; 6 M guanidine-HCl; 70% ethanol. 2 hours at room temperature: 0.1 M HCl, pH 1.0; 8 M urea, pH 10.5; 0.1 M Glycine-NaOH, pH 11.
Physical stability	Negligible volume variation due to changes in pH or ionic strength
pH stability ³	
Working range	3 to 9
Cleaning-in-place	2 to 10
Sanitization	Sanitize the column with 70% ethanol
Storage temperature	20% ethanol at 2°C to 8°C

¹ The dynamic capacity was calculated as the amount of human IgG bound to the resin before the flow through exceeded 1% of the absorbance of the in-going solution. The capacity was determined at a linear flow rate of 30 cm/h and with a sample concentration of 0.92 mg/mL. Please note that there might be considerable deviations in binding capacity for different immunoglobulins derived from the same species, even if they are of the same subclass.

² Linear flow rate = volumetric flow rate (cm³/h)/column cross-sectional area (cm²)

³ Working range: pH interval where the resin can be handled without significant change in function.

Cleaning-in-place: pH interval where the resin can be subjected to cleaning-in-place without significant change in function.

pH below 3 is sometimes required to elute strongly bound antibody species. However protein ligands may hydrolyze at very low pH.

Table 2. Relative binding strengths for protein A and protein G

Species	Subclass	Protein A binding ¹	Protein G binding
Human	IgA	variable	-
	IgD	-	-
	IgE	-	-
	IgG ₁	++++	++++
	IgG ₂	++++	++++
	IgG ₃	-	++++
	IgG ₄	++++	++++
	IgM	variable	-
Avian egg yolk	IgY	-	-
Cow		++	++++
Dog		++	+
Goat		-	++
Guinea pig	IgG ₁	++++	++
Hamster		+	++
Horse		++	++++
Koala		-	+
Llama		-	+
Monkey (rhesus)		++++	++++
Mouse ²	IgG ₁	+	++++
	IgG _{2a}	++++	++++
	IgG _{2b}	+++	+++
	IgG ₃	++	+++
	IgM	variable	-
Pig		+++	+++
Rabbit	No distinction	++++	+++
Rat ³	IgG ₁	-	+
	IgG _{2a}	-	++++
	IgG _{2b}	-	++
	IgG ₃	+	++
Sheep		+/-	++

¹ ++++ = strong binding; ++ = medium binding; - = weak or no binding

² IgG₃ from mouse binds more strongly to protein G than to protein A.

³ Note that IgG from rat binds to protein G coupled to Sepharose 4 Fast Flow.

3 Operation

Preparation of the chromatography resin

Protein G Sepharose 4 Fast Flow is supplied preswollen in 20% ethanol. Prepare a slurry by decanting the 20% ethanol solution and replace it with binding buffer in a ratio of 75% settled resin to 25% buffer. The binding buffer should not contain agents which significantly increase the viscosity. The column may be equilibrated with viscous buffers at reduced flow rate after packing is completed.

For batch procedures remove the ethanol by washing the resin on a medium porosity sintered glass funnel.

Packing of Protein G Sepharose 4 Fast Flow

Follow the instructions for how to pack Protein G Sepharose 4 Fast Flow into an empty column.

Step	Action
1	Equilibrate all material at the temperature at which the chromatography will be performed.
2	De-gas the resin slurry.
3	Eliminate air from the column dead spaces by flushing the end pieces with binding buffer. Make sure no air has been trapped under the column net. Close the column outlet with a few centimetres of binding buffer remaining in the column.
4	Pour the slurry into the column in one continuous motion. Pouring the slurry down a glass rod held against the wall of the column will minimise the introduction of air bubbles.
5	Immediately fill the remainder of the column with binding buffer. Mount the column top piece onto the column and connect the column to a pump.
6	Open the bottom outlet of the column and set the pump to run at the desired flow rate. ¹

Step	Action
------	--------

7	Maintain the packing flow rate for 3 bed volumes after a constant bed height is reached.
---	--

- ¹ Ideally, Fast Flow resins are packed at a constant pressure not exceeding 1 bar (0.1 MPa) in XK columns. If the packing equipment does not include a pressure gauge, use a packing flow rate of at least 400 cm/h (15 cm bed height, 25°C, low viscosity buffer).

If the recommended pressure or flow rate cannot be obtained, use the maximum flow rate the pump can deliver. This should also give a reasonably well-packed gel.

Note:

Do not exceed 75% of the packing flow rate in subsequent purification procedures.

Using an adapter

Follow the instructions to fit an adapter.

Step	Action
------	--------

- | | |
|---|--|
| 1 | After the resin has been packed as described above, close the column outlet and remove the top piece from the column. Carefully fill the rest of the column with distilled water to form an upward meniscus at the top. |
| 2 | Insert the adapter an angle into the column, ensuring that no air is trapped under the net. |
| 3 | Make all tubing connections at this stage. There must be a bubble-free liquid connection between the column and the pump and column. |
| 4 | Slide the plunger slowly down the column so that the air above the net and in the capillary tubings is displaced by buffer. |
| 5 | Lock the adapter in position, open the column outlet and start the flow of buffer. Pass buffer through the column at the packing flow rate until the resin bed is stable. Re-position the adapter on the resin surface as necessary. |

The column is now packed and equilibrated and ready for use.

Binding

IgG from most species binds Protein G Sepharose 4 Fast Flow at neutral pH and physiological ionic strength.

As a general method we recommend 20 mM sodium phosphate, pH 7.0 as binding buffer.

The binding capacity of Protein G Sepharose 4 Fast Flow depends on the source of the particular immunoglobulin, see [Table 3](#).

However, the total capacity depends upon several factors, such as the flow rate during sample application, the sample concentration and binding buffer. The table below shows the total capacity under defined conditions for IgG from some species.

Table 3. The total IgG capacity of Protein G Sepharose 4 Fast Flow, under defined conditions¹, for various species. (Work from GE)

Species	Total IgG capacity
Human	17
Rat	7
Sheep	18
Rabbit	19
Goat	19
Guinea-pig	17
Cow	23
Mouse ²	6

¹ The total capacity was determined with 1 mL drained resin packed in a 1 × 10 cm column at a linear flow rate of approx. 11 cm/h. Approximately 40 mg IgG was applied. Binding buffer used was 20 mM sodium phosphate, pH 7.0 and elution buffer used was 0.1 M glycine-HCl, pH 2.7.

² Extrapolated value from experiment carried out at 1/5th scale.

Note:

The binding capacity values listed above are typical for the given species. However, there might be considerable deviations in binding capacity for different immuno-globulins derived from the same species, even if they are of the same subclass.

Elution

To elute IgG from Protein G Sepharose 4 Fast Flow it is necessary to lower the pH to about 3.0 to 2.5 depending on the IgG.

As a general method, we recommend 0.1 M glycine buffer, pH 3.0 to 2.5 as elution buffer.

As a safety measure to preserve the activity of acid labile IgG's, it is recommended to add of 60 to 200 $\mu\text{L}/\text{mL}$ eluate of 1 M Tris-HCl, pH 9.0, to neutralize the eluted fractions.

Regeneration

After elution, the resin should immediately be washed with 2 to 3 bed volumes of elution buffer followed by re-equilibration with 2 to 3 bed volumes of binding buffer.

In some applications, substances like denatured proteins or lipids, do not elute in the regeneration procedure. These can be removed by cleaning-in-place procedures.

4 Cleaning-in-place and Sanitization

Cleaning-in-place (CIP)

Remove strongly bound hydrophobic proteins, lipoproteins and lipids by washing the column with a non-ionic detergent (e.g. Triton™ X-100), 0.1%, at 37°C, contact time one minute. Immediately re-equilibrate with at least 5 bed volumes of sterile filtered binding buffer.

Alternatively, wash the column with 70% ethanol and let stand for 12 hours. Re-equilibrate with at least 5 bed volumes of sterile binding buffer.

Sanitization

Sanitization reduces microbial contamination of the chromatography resin to a minimum.

Wash the column with a buffer containing 2% hibitane digluconate and 20% ethanol. Allow to stand for 6 hours.

Re-equilibrate the column with 3 to 5 bed volumes of sterile binding buffer.

Column performance is normally not significantly changed by the cleaning in place or sanitization procedures described above.

These recommended cleaning procedures can be performed directly on the packed column.

5 Immunoprecipitation

Immunoprecipitation is a highly specific and effective technique for analytical separations of target antigens from crude cell lysates.

Getting started

To obtain satisfactory results using immunoprecipitation, all procedures involved must be empirically optimized. For example, selecting cell lysis conditions is very critical and has to be optimized with regard to cell type and how the antigen is to be used. Whereas cells without cell walls (e.g., animal cells) are easily disrupted by treatment with mild detergent, other cells may need some type of mechanical shearing such as sonication or homogenization.

The parameters listed below (lysis buffers, incubation times, volumes, and concentrations) should therefore be regarded as guidelines for initial experiments.

For common problems and general tips, see instructions *Immunoprecipitation Starter Pack* (71-5017-54).

Preparing the resin

Protein G Sepharose 4 Fast Flow is supplied preswollen in 20% ethanol. Wash the resin three times with lysis buffer. Centrifuge at 12 000 g for 20 seconds between the washes and discard the supernatant. Prepare a 50% slurry by mixing equal volumes of resin and lysis buffer. Store at 4°C and mix well before use.

Cell lysis

Step	Action
1	<ul style="list-style-type: none">• Adherent cells: Remove all culture medium and wash twice with ice-cold PBS. Discard the supernatants and drain well.• Cells in suspension: Collect cells by centrifugation at 1000 g for 5 minutes and discard the culture medium supernatant . Resuspend the pellet in ice-cold PBS, centrifuge and discard the supernatant . Repeat the wash once.
2	<ul style="list-style-type: none">• Adherent cells: Place the tissue culture dish on ice. Add ice-cold lysis buffer¹ to a concentration of 10^6 to 10^7 cells/mL (1 mL to a cell culture plate, Ø 10 cm). Incubate on ice for 10 to 15 minutes with occasional rocking.• Cells in suspension: Suspend the washed pellet in ice-cold lysis buffer¹ at a concentration of 10^6 to 10^7 cells/mL (approximately 10 cell volumes lysis buffer). Incubate on ice for 10 to 15 minutes with gentle mixing.
3	Transfer the cells to a suitable homogenization tube.
4	Further disrupt the cells by sonication, homogenization or passage through a 21 Gauge needle. Keep the cells on an ice bath to prevent the temperature from rising.
5	Centrifuge at 12 000 g for 10 minutes at 4°C to remove particulate matter.
6	Transfer the lysate (the supernatant) to a fresh tube. Keep on ice.

¹ See Section *Buffers and solutions* for help when selecting lysis buffer.

Preclearing (optional)

Antibodies present in the cell lysate may also bind to the resin and thus interfere with subsequent analysis. In such a case preclearing may be desired.

Step	Action
1	Add 50 to 100 μ L Protein G Sepharose 4 Fast Flow suspension (50% slurry) to 1 mL cell lysate in an Eppendorf™ tube. Higher volume of resin might be necessary when working with serum samples due to the large amount of IgG present.
2	Gently mix for 1 hour at 4°C.
3	Centrifuge at 12 000 g for 20 seconds. Save the supernatant.

Couple antigen to antibody

Step	Action
1	Aliquot samples (500 μ L) in new Eppendorf tubes.
2	Add <ul style="list-style-type: none">• polyclonal serum (0.5 to 5 μL),• hybridoma tissue culture supernatant (5 to 100 μL),• ascites fluid (0.1 to 1 μL), or purified monoclonal or polyclonal antibodies (add the volume corresponding to 1 to 5 μg). <p>For controls, use non-immune antibodies that are as close to the specific antibody as possible (for example, polyclonal serum should be compared to normal serum from the same species).</p>
3	Gently mix for 1 hour at 4°C.

Precipitation of the immune complexes

Step	Action
1	Add 50 μ L Protein G Sepharose 4 Fast Flow suspension (50% slurry). Note: <i>It is possible to work with volumes down to 10 μL.</i>
2	Gently mix for 1 hour at 4°C.
3	Centrifuge at 12 000 g for 20 seconds and save the pellet.
4	Wash the pellet three times with 1 mL lysis buffer and once with wash buffer. Centrifuge at 12 000 g for 20 seconds between each wash and discard the supernatants. Note: <i>Be very careful when removing the supernatants to avoid loss of the beads!</i>

Dissociation and analysis

Step	Action
1	Suspend the final pellet in 30 μ L sample buffer.
2	Heat to 95°C for 3 minutes.
3	Centrifuge at 12 000 g for 20 seconds to remove the beads. Carefully remove the supernatant.
4	Add 1 μ L 0.1% bromphenol blue.
5	Analyze the supernatant by SDS-PAGE, followed by protein staining and/or immunoblotting for detection. Radiolabeled antigens are detected by autoradiography.

Buffers and solutions

Lysis buffers

Cell lysis must be harsh enough to release the target antigen, but mild enough to maintain its immunoreactivity. Selecting lysing conditions is therefore very critical and has to be individually optimized.

Some commonly used lysis buffers are listed below. NP-40 (IGEPAL CA-630) and RIPA buffer release most soluble cytoplasmic or nuclear proteins without releasing chromosomal DNA and are a good choice for initial experiments. Some parameters that affect the extraction of an antigen include salt concentration (0 to 1 M), non-ionic detergents (0.1 to 2%), ionic detergents (0.01 to 0.5%) and pH (6 to 9).

Name	Description	Stringency
Low salt	1% IGEPAL CA-630, 50 mM Tris, pH 8.0, 1 mM PMSF	+
Mammalian Protein Extraction Buffer	Tris-based buffer, 10 mM NaCl, detergent mixture (NP-40, Triton X-100, Tween™), pH 7.5	+
Yeast Protein Extraction Buffer	Tris-based buffer, 50 mM NaCl, detergent mixture (NP-40, Triton X-100, Tween), pH 7.5	+
NP-40 (IGEPAL CA-630)	150 mM NaCl, 1% IGEPAL CA-630, 50 mM Tris, pH 8.0, 1 mM PMSF	++
RIPA	150 mM NaCl, 1% IGEPAL CA-630, 0.5% sodium deoxycholate (DOC), 0.1% SDS, 50 mM Tris, pH 8.0, 1 mM PMSF	+++
High salt	500 mM NaCl, 1% IGEPAL CA-630, 50 mM Tris, pH 8.0, 1 mM PMSF	++++

Other buffers/solutions

Name	Description
PBS	1 mM KH_2PO_4 , 10 mM Na_2HPO_4 , 137 mM NaCl, 2.7 mM KCl, pH 7.4
Wash buffer	50 mM Tris, pH 8
Sample buffer (reducing)	1% SDS, 100 mM DTT, 50 mM, Tris, pH 7.5

6 Storage

For storage, keep the resin at 2°C to 8°C in a suitable bacteriostat, e.g., 20% ethanol. Protein G Sepharose 4 Fast Flow must not be frozen.

7 Ordering information

Product	Quantity	Code No.
Protein G Sepharose 4 Fast Flow ¹	5 mL	17061801
	25 mL	17061802
	200 mL	17061805
	1 L	17061806
	5 L	17061804

¹ Larger quantities are available. Contact GE for more information.

Related products	Quantity	Code No.
HiTrap™ Protein G HP	5 × 1 mL	17040401
	2 × 1 mL	17040403
	1 × 1 mL	29048581
	1 × 5 mL	17040501
	5 × 5 mL	17040503
MABTrap™ Kit	1 kit	17112801
Ab SpinTrap™	50 × 100 µL	28408347
Immunoprecipitation Starter Pack	2 × 2 mL	17600235
Protein G Sepharose 4 Fast Flow, 2 mL		
Protein G HP SpinTrap	16 × 100 µL	28903134
Protein G HP MultiTrap™	4 × 96-well plates	28903135
Ab Buffer Kit	1 kit	28903059
Protein G	5 mg	17061901
Protein A Mag Sepharose	500 µL	28944006
	4 × 500 µL	28951378
Protein G Mag Sepharose	500 µL	28944008
	4 × 500 µL	28951379

Product	Quantity	Code No
Antibody Purification Handbook	1	18103746
Solutions for antibody purification, Selections Guide	1	28935197
Affinity Chromatography Handbook	1	18102229
Convenient Protein Purification, HiTrap Column Guide	1	18112981
Prepacked chromatography columns for ÄKTA™ systems, Selection guide	1	28931778

For local office contact information, visit
www.gelifesciences.com/contact

GE Healthcare UK Limited
Amersham Place
Little Chalfont
Buckinghamshire, HP7 9NA
United Kingdom

www.gelifesciences.com/protein-purification

GE, GE monogram and ÄKTA, BioProcess, HiTrap, MAbTrap, MultiTrap, Sepharose and SpinTrap are trademarks of General Electric Company.

Eppendorf is a trademark of Eppendorf AG. Triton is a trademark of Union Carbide Chemicals and Plastic Company Inc. Tween is a trademark of Croda International PLC.

All other third party trademarks are the property of their respective owner.

© 2005-2016 General Electric Company
First published Mar. 2005

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Bio-Sciences AB
Björkgatan 30, 751 84 Uppsala, Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare Bio-Sciences Corp.
100 Results Way, Marlborough, MA 01752, USA

GE Healthcare Dharmacon, Inc.
2650 Crescent Dr., Lafayette, CO 80026, USA

HyClone Laboratories, Inc.
925 W 1800 S, Logan, UT 84321, USA

GE Healthcare Japan Corporation
Sanken Bldg. 3-25-1, Hyakunincho Shinjuku-ku, Tokyo 169-0073, Japan

