

# HiScreen™ Capto™ S ImpAct

## HiTrap™ Capto S ImpAct, 1 ml, 5 ml

Capto S ImpAct is a strong cation exchange BioProcess™ chromatography medium (resin) for intermediate purification and polishing of a wide range of biomolecules, especially monoclonal antibodies. The high binding capacity and the high flow base matrix in combination with the small particle size of Capto S ImpAct result in excellent pressure-flow properties as well as impressive resolution.

Capto S ImpAct provides:

- Efficient aggregate removal at high load of monoclonal antibodies
- High-resolution intermediate purification and polishing based on the well-established Capto platform with novel polymeric ligand
- Flexibility of design – large operational window of flow rates and bed heights
- High-throughput purifications easy to optimize and scale up
- Higher manufacturing productivity enables improved process economy

HiScreen Capto S ImpAct (4.7 ml) and HiTrap Capto S ImpAct (1 ml and 5 ml) are prepacked columns for small scale purifications, as well as optimization of methods and parameters, such as sample load and binding conditions. HiScreen Capto S ImpAct and HiTrap Capto S ImpAct columns provide fast, reproducible and easy separations in convenient formats. The columns are best utilized when they are connected to liquid chromatography systems such as ÄKTA™.



# Table of Contents

1	Product description .....	3
2	Optimization .....	7
3	Operation .....	8
4	Cleaning-in-place (CIP) .....	14
5	Scaling up .....	16
6	Adjusting pressure limits in chromatography system software .....	17
7	Storage .....	18
8	Troubleshooting .....	19
9	Ordering information .....	20

Please read these instructions carefully before using the products.

## **Intended use**

The products are intended for research use only, and shall not be used in any clinical or *in vitro* procedures for diagnostic purposes.


## **Safety**

For use and handling of the products in a safe way, please refer to the Safety Data Sheets.

# 1 Product description

## HiScreen column characteristics

HiScreen columns are made of biocompatible polypropylene that does not interact with biomolecules. The columns are delivered with stoppers at the inlet and outlet. The arrow on the column label indicates the column orientation and the recommended flow direction, see Figure 1.



**Fig 1.** HiScreen column

For easy scaling-up and when a higher bed height is required, two columns can easily be connected in series using a union to give 20 cm bed height (see Section *Scaling up*).

**Note:** *HiScreen columns cannot be opened or refilled*

**Note:** *Make sure that the connector is tight to prevent leakage.*

**Table 1.** Characteristics of HiScreen column

Column volume (CV)	4.7 ml
Column dimensions	0.77 × 10 cm
Column hardware pressure limit	8 bar (0.8 MPa)
Recommended flow rate	1.2 ml/min
Max flow rate	2.3 ml/min

**Note:** *The pressure over the packed bed varies depending on a range of parameters such as the characteristics of the chromatography medium, sample/liquid viscosity and the column tubing used.*

## HiTrap column characteristics

The columns are made of biocompatible polypropylene that does not interact with biomolecules.

The columns are delivered with a stopper at the inlet and a snap-off end at the outlet. Table 2 lists the characteristics of HiTrap columns.



Fig 2. HiTrap, 1 ml column.



Fig 3. HiTrap, 5 ml column

**Note:** *HiTrap columns cannot be opened or refilled.*

**Note:** *Make sure that the connector is tight to prevent leakage.*


**Table 2.** Characteristics of HiTrap columns

Column volume (CV)	1 ml	5 ml
Column dimensions	0.7 × 2.5 cm	1.6 × 2.5 cm
Column hardware pressure limit <sup>1</sup>	5 bar (0.5 MPa)	5 bar (0.5 MPa)
Recommended flow rate	1 ml/min	5 ml/min
Max flow rate	4 ml/min	20 ml/min

<sup>1</sup> The pressure over the packed bed varies depending on a range of parameters such as the characteristics of the chromatography medium and the column tubing used.

## Properties of Capto S ImpAct

The structure of the well established S (sulfonate group) ligand used for Capto S ImpAct is shown in Figure 4.



**Fig 4.** Sketch showing the strong ion exchange group of Capto S ImpAct

The small bead size of Capto S ImpAct allows for high resolution purifications. The combination of the well established S ligand with a small high flow agarose bead makes the media ideal for intermediate purification and polishing. Further characteristics of the media are found in Table 3.

### Polishing

Capto S ImpAct is designed for post protein A polishing of monoclonal antibodies. Removal of contaminants such as aggregates, host cell proteins, and leached protein A in either bind-elute or flow-through mode.

### High flow velocities

In spite of the small bead size of Capto S ImpAct, the rigid high flow agarose matrix allows high flow velocities.

## Characteristics

The highly cross-linked agarose base matrix gives the medium high chemical and physical stability. Characteristics such as capacity, elution behavior and pressure/flow properties are unaffected by the solutions commonly used in process chromatography and cleaning procedures (see table below).

**Table 3.** Characteristics of Capto S ImpAct

<b>Matrix</b>	High flow agarose
<b>Functional group</b>	SO <sub>3</sub> <sup>-</sup>
<b>Total ionic capacity</b>	37 to 63 µmol (H <sup>+</sup> ) /ml medium
<b>Average particle size (d<sub>50v</sub>)<sup>1</sup></b>	50 µm
<b>Maximum operational flow velocity<sup>2</sup></b>	At least 220 cm/h in a 1 m diameter column with bed height 20 cm at 20°C; measured using process buffers with the same viscosity as water at < 3 bar (0.3 MPa)
<b>Binding capacity<sup>3</sup> (mg/ml medium)</b>	> 90 mg lysozyme > 85 mg BSA > 100 mg IgG
<b>pH stability<sup>4</sup></b>	
Working	4 to 12
Cleaning-in-place	3 to 14
<b>Working temperature</b>	4°C to 30°C
<b>Chemical stability</b>	All commonly used aqueous buffers, 1 M sodium hydroxide <sup>5</sup> , 8 M urea, 6 M guanidine hydrochloride, 30% isopropanol and 70% ethanol
<b>Avoid</b>	Oxidizing agents, cationic detergents
<b>Storage</b>	0.2 M sodium acetate in 20% ethanol

<sup>1</sup> d<sub>50v</sub> is the median particle size of the cumulative volume distribution.

<sup>2</sup> Specific flow values for each column type are shown in Table 1 (HiScreen) and Table 2 (HiTrap).

<sup>3</sup> Dynamic binding capacity at 10% breakthrough measured at a residence time of 5.4 minutes (220 cm/h) in a Tricorn™ 5/50 column with 5 cm bed height. 25 mM sodium phosphate, pH 7.2 (lysozyme), 50 mM sodium acetate, pH 5.0 (BSA) and 50 mM sodium acetate, pH 5.5 (IgG)

<sup>4</sup> Working range: pH interval where the medium can be operated without significant change in function. Cleaning-in-place: pH stability where the medium can be subjected to cleaning-in-place without significant change in function.

<sup>5</sup> No significant change in ionic capacity and carbon content after 1 week storage in 1 M NaOH at 40°C.

## 2 Optimization

### Aim

The aim of designing and optimizing an ion exchange separation process of biomolecules is to identify conditions that promote binding of the highest amount of target molecule, in the shortest possible time with highest possible product recovery and purity. For optimization of binding conditions, pH and conductivity (salt concentration) should be screened.

### Workflow

The recommended workflow is described in the figure below. It starts with screening of conditions in high throughput formats, followed by optimization, preferably by applying Design of Experiments (DoE), in small columns and finally scale-up to large columns.



Fig 5. Recommended workflow during process development.

### Prepacked chromatography media

Prefilled PreDicator™ plates or PreDicator RoboColumn™ are recommended for condition screening while prepacked HiScreen™ columns are recommended for optimization. Optimization in columns is preferentially done by using Design of Experiments (DoE). DoE can also be used to evaluate plate data, but often raw data analysis of plate data is enough to find conditions that can be used during optimization in packed bed. For experimental setup, data handling, analysis and evaluation of raw data from plates, the Assist software from GE Healthcare will be useful.

# 3 Operation

## Prepare buffers

The buffer species and buffer concentration are important for reproducible and robust methods. The table below shows suitable buffers for ion exchange chromatography and suggested starting concentrations.

Users of ÄKTA systems with BufferPrep or BufferPro functionality can select from a range of buffer recipes to conveniently screen media over a range of pH values and conductivities.

**Table 4.** Recommended buffers for ion exchange chromatography

pH interval	Buffer <sup>1, 2</sup>	Concentration <sup>3</sup>
4 to 5.5	Acetate	20 to 100 mM
3 to 6.5	Citrate	20 to 200 mM
5.5 to 6.5	Bis-Tris	20 to 50 mM
6 to 7.5	Phosphate	50 to 200 mM
7.5 to 9	Tris	20 to 50 mM
8.8 to 10.6	Glycin-NaOH	20 to 100 mM

<sup>1</sup> The choice of buffer systems and salts may influence both yield and contaminant clearance.

<sup>2</sup> Buffers in the interval 5.5 to 8 will normally be most efficient for contaminant removal.

<sup>3</sup> Conductivity can further be adjusted by addition of salt.

## Prepare the sample

Step	Action
1	Adjust the sample to the composition of the start buffer, using one of these methods: <ul style="list-style-type: none"><li>• Dilute the sample with start buffer.</li><li>• Exchange buffer using a HiPrep™ 26/10 Desalting, HiTrap™ Desalting or PD-10 Desalting column (see table below).</li></ul>
2	Filter the sample through a 0.45 µm filter or centrifuge immediately before loading it to the column. This prevents clogging and increases the life time of the column when loading large sample volumes.



**Note:** *For larger volumes of feed, sample preparation is preferably performed by diafiltration or directly by adjustment of pH and conductivity.*

**Table 5.** Prepacked columns for desalting

Column	Code No.	Loading volume	Elution volume	Comments	Application
HiPrep 26/10 Desalting	17-5087-01	2.5 to 15 ml	7.5 to 20 ml	Prepacked with Sephadex™ G-25 Fine. Requires a laboratory pump or a chromatography system to run.	For desalting and buffer exchange of protein extracts ( $M_r > 5000$ ).
HiTrap Desalting	17-1408-01	0.25 to 1.5 ml	1.0 to 2.0 ml	Prepacked with Sephadex G-25 Superfine. Requires a syringe or pump to run.	
PD-10 Desalting	17-0851-01	1.0 to 2.5 ml <sup>1</sup> 1.75 to 2.5 ml <sup>2</sup>	3.5 ml <sup>1</sup> up to 2.5 ml <sup>2</sup>	Prepacked with Sephadex G-25 Medium.	For desalting, buffer exchange, and cleanup of proteins and other large biomolecules ( $M_r > 5000$ ).
PD MiniTrap™ G-25	28-9180-07	0.1 to 0.5 ml <sup>1</sup> 0.2 to 0.5 ml <sup>2</sup>	1.0 ml <sup>1</sup> up to 0.5 ml <sup>2</sup>	Runs by gravity flow or centrifugation	
PD MidiTrap™ G-25	28-9180-08	0.5 to 1.0 ml <sup>1</sup> 0.75 to 1.0 ml <sup>2</sup>	1.5 ml <sup>1</sup> up to 1.0 ml <sup>2</sup>		

<sup>1</sup> Volumes with gravity elution

<sup>2</sup> Volumes with centrifugation

## Recommended flow rates

The pressure-flow specification for the Capto S ImpAct medium in large columns with negligible wall support is >220 cm/h in a 20 cm bed height with a buffer viscosity comparable to the one of water. Due to pronounced wall support for the prepacked formats, higher flow rates are possible to use in equilibrium and wash phases. The residence time used for sample application, elution and CIP should however match those appropriate for the final intended scale of the process. Recommended flow rates for different operations are found in the tables below.

**Table 6.** Recommended flow rates using HiScreen Capto S ImpAct

Type of operation	Flow rate (ml/min)	Flow velocity (cm/h)	Residence time (min)
Equilibration <sup>1</sup>	≤ 2.3	≤ 300	≥ 2
Wash <sup>1</sup>	≤ 2.3	≤ 300	≥ 2
Sample load	0.78 to 1.55	100 to 200	3 to 6
Cleaning-in-place <sup>2</sup>	≤ 0.78	≤ 100	≥ 6

**Table 7.** Recommended flow rates using HiTrap Capto S ImpAct, 1 ml

Type of operation	Flow rate (ml/min)	Flow velocity (cm/h)	Residence time (min)
Equilibration <sup>1</sup>	≤ 4.0	≤ 600	≥ 0.25
Wash <sup>1</sup>	≤ 4.0	≤ 600	≥ 0.25
Sample load	0.16 to 0.32	25 to 50	3 to 6
Cleaning-in-place <sup>2</sup>	≤ 0.16	≤ 25	≥ 6

**Table 8.** Recommended flow rates using HiTrap Capto S ImpAct, 5 ml

Type of operation	Flow rate (ml/min)	Flow velocity (cm/h)	Residence time (min)
Equilibration <sup>1</sup>	≤ 20	≤ 600	≥ 0.25
Wash <sup>1</sup>	≤ 20	≤ 600	≥ 0.25
Sample load	0.84 to 1.68	25 to 50	3 to 6
Cleaning-in-place <sup>2</sup>	≤ 0.84	≤ 25	≥ 6

<sup>1</sup> Stated flow rates are for buffers with the same viscosity as water at 20°C. For solutions with higher viscosities, e.g., 20% ethanol, lower flow rates should be used.

<sup>2</sup> The Cleaning-in-place should be done in at least 3 column volumes with a total contact time of at least 15 min. See also *Chapter 4 Cleaning-in-place (CIP)*, on page 14.

## Purification

Collect fractions throughout the separation.

**Flow rate:** See the tables 6, 7 and 8 for recommended flow rates for different operations and formats.

**Column tubing:** Choose the optimal tubing kit for the column and the application you intend to run. (i.d.: 0.25, 0.50 or 0.75 [mm]). A tubing with wider inner diameter gives broader peaks whereas a tubing with a smaller inner diameter gives higher back pressure.

### Column and sample preparation

Step	Action
------	--------

- |   |   |
|---|---|
| 1 | Before connecting the column to a system: <ul style="list-style-type: none"><li>• Remove the two stoppers from the HiScreen column.</li><li>• Remove the stopper in the top and the snap-off in the bottom of the HiTrap column.</li></ul>  |
| 2 | Connect the column to the system with a drop-to-drop connection to avoid introducing air into the column.<br><b>Note:</b><br><i>To prevent leakage, make sure that the connectors are tight. Use fingertight 1/16" connector (28-4010-81).</i>  |
| 3 | Wash with 1 column volume (CV) of distilled water. This step removes the ethanol and avoids precipitation of buffer salts upon exposure to ethanol. Since 20% ethanol has higher viscosity than water, the recommended flow rate for this step is: <ul style="list-style-type: none"><li>• 1.2 ml/min for HiScreen (150 cm/h)</li><li>• 1.0 ml/min for HiTrap, 1 ml (150 cm/h)</li><li>• 5.0 ml/min for HiTrap, 5 ml (150 cm/h)</li></ul> <b>Note:</b><br><i>This step can be omitted if precipitation is not likely to be a problem.</i> |
| 4 | Equilibrate the column with at least 5 CV start buffer or until the UV baseline, eluent pH and conductivity are stable.   |

Step	Action
5	Adjust the sample to the chosen starting pH and conductivity, see table 5.
6	Load sample onto the column. Collect fractions for analysis of yield and purity. <b>Note:</b> <i>Target protein is bound.</i>
7	Wash out unbound material with start buffer.
8	Elute bound material and collect elution fractions for analysis of yield and purity.
9	Wash with 1 M NaCl (100% elution buffer) to elute any remaining ionically bound material.
10	If required, perform a CIP to elute and clean the column, see Section <i>Cleaning-in-place (CIP)</i> .
11	Equilibrate with 5-10 CV starting buffer or until the UV baseline, pH, and conductivity reach the required values. The equilibration step can be shortened by first washing with a high concentration buffer to obtain approximately the desired pH value and then washing with starting buffer until the conductivity and pH values are stable. <b>Note:</b> <i>Do not exceed the maximum recommended flow velocity or the column hardware pressure limit (see Tables 1 and 2).</i>

---

## 4 Cleaning-in-place (CIP)

### General description

Correct preparation of samples and buffers should maintain columns in good condition. However, reduced performance, increased back pressure or blockage indicates that the column medium needs cleaning.

CIP removes very tightly bound, precipitated or denatured substances from the medium. If such contaminants are allowed to accumulate, they may affect the chromatographic properties of the prepacked column, reduce the capacity of the medium and, potentially, come off in subsequent runs. If the fouling is severe, it may block the column, increase back pressure and affect the flow properties.

CIP should be performed regularly to prevent the build-up of contaminants and to maintain the capacity, flow properties and general performance of prepacked columns. CIP is usually performed after the elution.

It is recommended to perform a CIP:

- When an increase in back pressure is seen.
- If reduced column performance is observed.
- Before first time use or after long term storage.
- Between runs when the same column is used for purification of different proteins to prevent possible cross-contamination.
- After every run with real feed.

## CIP protocol

The nature of the sample will ultimately determine the final CIP protocol so the CIP procedure below may require optimization. NaOH concentration, contact time and frequency are typically the main parameters to vary during the optimization of the CIP.

The CIP procedure below removes common contaminants.

### Flow rate

HiScreen	0.78 ml/min (100 cm/h)
HiTrap, 1 ml	0.16 ml/min (25 cm/h)
HiTrap, 5 ml	0.84 ml/min (25 cm/h)

Step	Action
1	Wash with at least 2 column volumes (CV) of 2 M NaCl.
2	Wash with at least 3 CV 1 M NaOH with at least 15 min contact time.
3	Wash with at least 2 CV 2 M NaCl.
4	Wash with at least 2 CV distilled water.
5	Wash with 5 CV start buffer or until eluent pH and conductivity have reached the required values

For further information about CIP, see instruction for Capto S ImpAct, 29-0925-01.

## 5 Scaling up

After optimizing the method at laboratory scale, the process is ready for scaling up. For quick small scale-up of purification, either two HiTrap columns (one in the top of the other) or two HiScreen columns can conveniently be connected in series with a union (18-1120-93) to give increased bed height. Note that the back pressure will increase. This can easily be addressed by lowering the flow rate.

Factors such as clearance of critical impurities may change when column bed height is modified and should be validated using the final bed height.

Scaling up is otherwise typically performed by keeping bed height and linear flow velocity (cm/h) constant while increasing bed diameter and volumetric flow rate (ml/min or l/h).

Bulk medium is available for further scaling up, see *Ordering information*.

A general description of the scaling up procedure is described below.

<b>Step</b>	<b>Action</b>
<b>1</b>	Select bed volume according to required sample load. Keep sample concentration constant.
<b>2</b>	Select column diameter to obtain the desired bed height. The excellent rigidity of the high flow base matrix allows for flexibility in choice of bed heights.
<b>3</b>	The larger equipment used when scaling up may cause some deviations from the method optimized at small scale. In such cases, check the buffer delivery and monitoring systems for time delays or volume changes.




## 6 Adjusting pressure limits in chromatography system software

Pressure generated by the flow through a column affects the packed bed and the column hardware, see Figure below. Increased pressure is generated when running/using one or a combination of the following conditions:

- High flow rates
- Buffers or sample with high viscosity
- Low temperature
- A flow restrictor

**Note:** *Exceeding the flow limit (see Table 1) may damage the column.*



**Fig 6.** Pre-column and post-column measurements.

### ÄKTA avant

The system will automatically monitor the pressures (pre-column pressure and pressure over the packed bed,  $\Delta p$ ). The pre-column pressure limit is the column hardware pressure limit (see Tables 1 and 2).

The maximum pressure the packed bed can withstand depends on medium characteristics and sample/liquid viscosity. The measured value also depends on the tubing used to connect the column to the instrument.

## ÄKTAexplorer, ÄKTApurifier, ÄKTAFPLC and other systems with pressure sensor in the pump

To obtain optimal functionality, the pressure limit in the software may be adjusted according to the following procedure:

Step	Action
1	Replace the column with a piece of tubing. Run the pump at the maximum intended flow rate. Note the pressure as <i>total system pressure</i> , P1.
2	Disconnect the tubing and run the pump at the same flow rate used in step 1. Note that there will be a drip from the column valve. Note this pressure as P2.
3	Calculate the new pressure limit as a sum of P2 and the column hardware pressure limit (see Tables 1 and 2). Replace the pressure limit in the software with the calculated value. The actual pressure over the packed bed ( $\Delta p$ ) will during run be equal to actual measured pressure - <i>total system pressure</i> (P1).

**Note:** Repeat the procedure each time the parameters are changed.

## 7 Storage

Wash with 2 column volumes (CV) of distilled water followed by 2 CV 20% ethanol, 0.2 M sodium acetate.

Store at 4°C to 30°C. Do not freeze.

Make sure that the column is tightly sealed to avoid drying out.

## 8 Troubleshooting

<b>Problem</b>	<b>Possible cause/corrective action</b>
High back pressure during the run	The column is clogged. <i>Clean the column, see Section Cleaning-in-place (CIP).</i> High viscosity of solutions. <i>Use lower flow rate.</i>
Unstable pressure curve during sample loading	Air bubbles trapped in the sample pump. <i>If possible, degas the sample using a vacuum degasser.</i>
Gradual decrease in yield	Insufficient elution and CIP. <i>Optimize the elution conditions, the CIP protocol and/or perform CIP more frequently.</i>
Gradual increase in CIP peaks	Sub-optimal elution conditions and CIP. <i>Optimize the elution conditions, the CIP protocol and/or perform CIP more frequently.</i>
Reduced column performance despite optimized elution and CIP	<i>Change to a new column. The longevity of the column depends mainly on the sample and sample preparation.</i>

## 9 Ordering information

<b>Product</b>	<b>Quantity</b>	<b>Code No.</b>
HiScreen Capto S ImpAct	1 × 4.7 ml	17-3717-47
HiTrap Capto S ImpAct	5 × 1 ml	17-3717-51
HiTrap Capto S ImpAct	5 × 5 ml	17-3717-55

<b>Related products</b>	<b>Quantity</b>	<b>Code No.</b>
Capto S ImpAct	25 ml	17-3717-01
	100 ml	17-3717-02
	1 l	17-3717-03
	5 l	17-3717-04
	10 l	17-3717-05
HiTrap Desalting	1 × 5 ml	29-0486-84
	5 × 5 ml	17-1408-01
HiPrep 26/10 Desalting	1 × 53 ml	17-5087-01
	4 × 53 ml	17-5087-02

<b>Accessories HiScreen</b>	<b>Quantity</b>	<b>Code No.</b>
HiTrap/HiPrep, 1/16" male connector for ÄKTA <i>(For connection of columns with 1/16" fittings to ÄKTA)</i>	8	28-4010-81
Union 1/16" male/1/16" male with 0.5 mm i.d. <i>(For connecting two columns with 1/16" fittings in series)</i>	2	18-1120-93
Fingertight stop plug, 1/16" <sup>1</sup> <i>(For sealing a HiScreen column)</i>	5	11-0003-55

<sup>1</sup> One fingertight stop plug is connected to the inlet and the outlet of each HiScreen column at delivery.

<b>Accessories HiTrap</b>	<b>Quantity</b>	<b>Code No.</b>
1/16" male/luer female <i>(For connection of syringe to top of HiTrap column)</i>	2	18-1112-51
Tubing connector flangeless/M6 female <i>(For connection of tubing to bottom of HiTrap column)</i>	2	18-1003-68
Tubing connector flangeless/M6 male <i>(For connection of tubing to top of HiTrap column)</i>	2	18-1017-98
Union 1/16" female/M6 male <i>(For connection to original FPLC System through bottom of HiTrap column)</i>	6	18-1112-57
Union M6 female/ 1/16" male <i>(For connection to original FPLC System through top of HiTrap column)</i>	5	18-3858-01
Union luerlock female/M6 female	2	18-1027-12
HiTrap/HiPrep, 1/16" male connector for ÄKTA design	8	28-4010-81
Stop plug female, 1/16" <i>(For sealing bottom of HiTrap column)</i>	5	11-0004-64
Fingertight stop plug, 1/16"	5	11-0003-55

<b>Related literature</b>	<b>Code No.</b>
Data file, Capto S ImpAct	29-0670-18
Handbook: Ion Exchange Chromatography & Chromatofocusing: Principles and Methods	11-0004-21
Application note: Column efficiency testing	28-9372-07

Page intentionally left blank

Page intentionally left blank

For local office contact information, visit  
[www.gelifesciences.com/contact](http://www.gelifesciences.com/contact)

GE Healthcare Bio-Sciences AB  
Björkgatan 30  
751 84 Uppsala  
Sweden

[www.gelifesciences.com/protein-purification](http://www.gelifesciences.com/protein-purification)

GE and GE monogram are trademarks of General Electric Company.

ÄKTA, AxiChrom, BioProcess, Capto, HiScreen, HiTrap and PreDictor are trademarks of GE Healthcare companies.

© 2014 General Electric Company – All rights reserved.  
First published May 2014

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Europe GmbH  
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Limited  
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp.  
800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Japan Corporation  
Sanken Bldg. 3-25-1, Hyakunincho Shinjuku-ku, Tokyo 169-0073, Japan

