



Affinity chromatography

Capto™ Lentil Lectin

Capto Lentil Lectin is an affinity chromatography medium (resin) for purification of glycoproteins and other molecules containing carbohydrates such as α -D-mannose and α -D-glucose or sterically related residues. The product is based on a rigid high-flow agarose base matrix that enables high flow rates. Capto Lentil Lectin medium is available in bulk as well as in various formats for high-throughput process development (Fig 1).

Key features of Capto Lentil Lectin include:

- Group-specific adsorbent for carbohydrate-containing molecules
- High productivity and cost-efficiency in downstream operations
- Animal component-free production
- Security of supply and regulatory support

Characteristics

Lectins are proteins that interact specifically and reversibly with certain carbohydrate residues. Immobilized lectins are valuable tools for isolating and separating glycoproteins, glycolipids, polysaccharides, subcellular particles, and for purifying detergent-solubilized cell membrane components. Lectins are also useful for assessing changes in levels or composition of surface glycoproteins during cell development and in malignant or virally transformed cell variants.

Capto Lentil Lectin medium is designed by coupling lentil lectin to an N-hydroxysuccinimide (NHS)-activated high-flow agarose base matrix (Fig 2). Lentil lectin is a metalloprotein isolated from *Lens culinaris* (lentil) seeds and binds to molecules containing α -D-mannose, α -D-glucose, or sterically related residues. To maintain the binding characteristics of Capto Lentil Lectin, the presence of both Mn^{2+} and Ca^{2+} is essential. These ions are present in large excess in the recommended storage solution for the medium. The lectin-metal ion complex remains active and is stable at neutral pH, even in the absence of free metal ions. However, to preserve the binding activity of the medium at pH < 5, excess Mn^{2+} and Ca^{2+} (1 mM) should be present in the buffer solutions. Main characteristics of Capto Lentil Lectin are summarized in Table 1.



Fig 1. Capto Lentil Lectin medium is produced in an animal component-free environment and is covered by our security of supply program.




Fig 2. Schematic structure of Capto Lentil Lectin medium.

Table 1. Main characteristics of Capto Lentil Lectin medium

Matrix	Highly cross-linked agarose
Average particle size (d_{50v})*	75 μ m
Ligand	Lentil (<i>Lens culinaris</i>) lectin
Ligand density	3 g/L
Coupling chemistry	NHS
Binding capacity†	~ 15 mg porcine thyroglobulin/mL medium
pH stability‡	3 to 10
Chemical stability	Stable in all commonly used aqueous buffers.
Avoid	Chelating agents such as EDTA, 8 M urea, or solutions having a pH < 3 should be avoided, as these conditions result in removal of manganese from the lectin, leading to loss of activity of the medium.
Flow velocity§	100 to 300 cm/h
Storage	2°C to 8°C in 20% ethanol containing 150 mM NaCl, 1 mM $CaCl_2$ and 1 mM $MnCl_2$

* d_{50v} is the median particle size of the cumulative volume distribution.

† Static binding capacity determined in 0.1 M phosphate buffer, pH 7.0

‡ pH interval where the medium can be operated and be subjected to cleaning in place without significant change in function.

§ Maximum recommended operating velocity for Capto Lentil Lectin is 300 cm/h at 20 cm bed height, determined with water at 20°C. Operating pressure is less than 2 bar.

Operation

Binding

Binding of glycoproteins and other carbohydrate-containing molecules to Capto Lentil Lectin medium occurs at a neutral pH in the presence of both Mn^{2+} and Ca^{2+} . These ions are present in excess in the solution in which the medium is supplied. The protein-metal ion complex remains active and is stable at neutral pH even in the absence of the free metal ions. However, to preserve the binding activity at $pH < 5$, excess Mn^{2+} and Ca^{2+} (1 mM) is required. Recommended binding buffer is 20 mM Tris-HCl, pH 7.4 containing up to 0.5 M NaCl to avoid non-specific ionic interactions.

Elution

Elution of bound substances can be achieved using an increasing gradient (continuous or step) of α -D-methylmannoside or α -D-methylglucoside. These carbohydrates act as strong eluents and many substances elute at 0.1 to 0.2 M. Higher concentrations might be required for more tightly bound substances. Glucose and mannose may also be used, but are weaker eluents. Strongly bound substances can also be eluted using low pH (within operating range) or with a 0.1 M borate buffer, pH 6.5. Elution of strongly bound substances can be facilitated by including 1% deoxycholate or other detergent in the elution buffer.

Regeneration

Capto Lentil Lectin can be regenerated by washing the medium with two to three bed volumes of a buffer solution containing 0.5 M NaCl, alternately with high pH (8.5) and low pH (5.5) between wash cycles. These cycles should be repeated three times followed by re-equilibration with three to five bed volumes of binding buffer. All strongly bound substances might not elute during the regeneration procedure. In such cases, a borate buffer containing 0.1% non-ionic detergent could be used at a low flow rate. A 20% ethanol wash or a gradient wash with up to 50% ethylene glycol may also be used to elute strongly bound substances.

As an alternative regeneration method, the medium can be washed with a detergent solution (e.g., 0.1% Triton™ X-100) at 37°C for 1 min. Re-equilibrate with at least five bed volumes of binding buffer after regeneration.

Chemical stability

The chemical stability of Capto Lentil Lectin was determined by a total organic carbon/nitrogen (TOC/TN) leakage analysis after storage in buffers of various pH for one week at 40°C (Fig 3). The results show that Capto Lentil Lectin is stable between pH 3 and 10. At a $pH < 3$ or $pH > 10$, the leakage of both carbon nitrogen increases.




Fig 3. Relative loss of carbon and nitrogen from Capto Lentil Lectin medium when stored in different pH for one week at 40°C.

Storage

Capto Lentil Lectin is supplied preswollen as a suspension in 20% ethanol containing 150 mM NaCl, 1 mM $CaCl_2$, and 1 mM $MnCl_2$ (storage solution). Recommended storage is at 2°C to 8°C in storage solution.

Ordering information

Product	Quantity	Product code
Capto Lentil Lectin	25 mL	17548901
Capto Lentil Lectin	100 mL	17548902
Capto Lentil Lectin	1 L	17548903
Capto Lentil Lectin	5 L	17548904
HiTrap™ Capto Lentil Lectin	5 × 1 mL	17548911
HiTrap Capto Lentil Lectin	1 × 5 mL	17548912
HiScreen™ Capto Lentil Lectin	2 × 4.7 mL	29157958

Related literature	Product code
HiScreen prepacked columns, data file	28930581
Affinity Chromatography Handbook, Principle and Methods	18102229

www.gelifesciences.com/bioprocess

GE and GE monogram are trademarks of General Electric Company. BioProcess, Capto, HiScreen, and HiTrap are trademarks of General Electric Company or one of its subsidiaries. Triton is a trademark of Union Carbide Chemicals and Plastic Company Inc. All other third party trademarks are the property of their respective owner.
© 2015 General Electric Company—All rights reserved. First published Jun. 2015
All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.
GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK
GE Healthcare Europe, GmbH, Munzinger Strasse 5, D-79111 Freiburg, Germany
GE Healthcare Bio-Sciences Corp., 800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA
GE Healthcare Japan Corporation, Sanken Bldg., 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan
For local office contact information, visit www.gelifesciences.com/contact