

Ready-To-Go™ RT-PCR Beads

Description

Ready-To-Go RT-PCR Beads are stable at room temperature and designed for performing single-tube RT-PCR. The reagents are optimized for full-length cDNA synthesis to >7.5 kb and optimal sensitivity from PCR.

Applications

Ready-To-Go RT-PCR Beads can be used with either total RNA or mRNA isolated using various methods (Fig 1). RT-PCR can be applied to a variety of analyses, from detecting RNA in a species to determining relative levels of a specific RNA (Fig 2). Gene-specific PCR primers can be added at the same time as the cDNA synthesis primer to perform one-step RT-PCR, or PCR primers can be added after first-strand synthesis for two-step reactions. The beads are available in either 0.5 ml or 0.2 ml tubes—enough for 100 or 96 reactions, respectively.

Ready-To-Go RT-PCR Beads offer:

- **Convenience:** All the components needed for RT-PCR are included in a reaction bead. Simply add water, template RNA, and primers.
- **Reproducibility:** Beads are manufactured under highly controlled conditions, ensuring reproducible results (Fig 3).
- **Reduced risk of contamination:** Individual predisposed reactions minimize sample handling and pipetting steps, thus reducing the risk of RNA degradation, contamination, and pipetting errors.
- **Reliability:** Each batch of beads is function-tested to ensure quality and consistency of results.
- **Ambient temperature stability:** Ready-To-Go RT-PCR Beads can be shipped and stored at room temperature.




Fig 1. RT-PCR detection of specific mRNAs from a range of sources. Source organisms and RNAs are indicated above the lanes. The target genes for amplification and the amounts of RNA used are indicated in the above table.



Fig 2. Relative titer of hemolymph juvenile hormone binding protein mRNA (6). Duplicate one-step RT-PCRs were run using 100 ng and 2 ng of total RNA for each time point. Five µg of pd(N)₆ and 10 pmol of the gene-specific primers were used for all reactions. Total RNA was isolated from a pool of 10 to 20 eggs from the black larval mutant of *Manduca sexta*.






Fig 3. Reproducibility of RT-PCR using Ready-To-Go RT-PCR Beads. A 450 bp fragment from the mouse mitotic checkpoint regulator (3) gene was amplified by two researchers (A and B) using three batches of beads.

Components

RT-PCR Beads: Room-temperature stable beads containing buffer, dATP, dCTP, dGTP, dTTP, approximately 2.0 units of *Taq* DNA Polymerase, M-MuLV reverse transcriptase, RNase inhibitor, and stabilizers (including RNase/ DNase-Free BSA). When a bead is reconstituted to 50 μ l, the concentration of each dNTP is 200 μ M in 10 mM Tris-HCl (pH 9.0), 60 mM KCl, and 1.5 mM MgCl₂.

Control Mix Beads: A bead containing rabbit globin mRNA (1 ng) and 8 pmol each of 5'-specific globin primer (5'-d[ACACTCTGGTCCATCCGACTGAG]-3') and 3'-specific globin primer (5'-d[GCCACTCACTCAGACTTTATTCAA]-3').

pd(N)₆: Lyophilized; 275 mg.

pd(T)₁₂₋₁₈: Lyophilized; 55 mg.

Quality Control

Each batch of Ready-To-Go RT-PCR Beads is tested to ensure its ability to generate an RT-PCR product using a Control Mix Bead and an RT-PCR bead. In addition, the beads are tested using *Drosophila melanogaster* mRNA (7.6 kb) and primers specific for a 425 bp sequence from the *fsh* gene.

Storage

At ambient temperature in pouch with desiccant.

References

1. Chomczynski, P. and Sacchi, N. Single-step method of RNA isolation by acid guanidinium thiocyanate-phenol-chloroform extraction. *Anal. Biochem.* **162**, 156–159 (1987).
2. Stayton M. M. *et al.*, Characterization of a full-length petunia cDNA encoding a polypeptide of the light-harvesting complex associated with photosystem I. *Plant Mol. Biol.* **10**, 127–137 (1987).
3. Starborg M. *et al.*, A novel murine gene encoding a 216-kDa protein is related to a mitotic checkpoint regulator previously identified in *Aspergillus nidulans*. *J. Biol. Chem.* **269**, 24133–24137 (1994).
4. Haynes, S. R. *et al.*, The *Drosophila* fsh locus, a maternal effect homeotic gene, encodes apparent membrane proteins. *Dev. Biol.* **134**, 246–257 (1989).
5. Lamb, P. and Crawford, Characterization of the human p53 gene. L., *Mol. and Cell. Biol.* **6**, 1379–1385 (1986).
6. Goodman, W. G. and Chang, E.S. Juvenile hormone cellular and hemolymph binding protein, in *Comprehensive Insect Physiology, Biochemistry, and Pharmacology: Endocrinology I*. (Kerkut, G. A. and Gilbert, L. I., eds.) pp 491, Pergamon Press, Oxford, England (1985).

Ordering information

Product	Code number
Ready-To-Go RT-PCR Beads (0.5 ml tubes) 100 reactions	27-9266-01
Ready-To-Go RT-PCR Beads (0.2 ml tubes) 96 reactions	27-9267-01
Ready-To-Go RT-PCR Beads (0.2 ml hinged tube with cap) 96 reactions	27-9259-01

Related products

First-strand synthesis

First-Strand cDNA Synthesis Kit	27-9261-01
Ready-To-Go You-Prime First-Strand Beads	27-9264-01
Ready-To-Go T-Primed First-Strand Kit	27-9263-01

PCR

Product	Code number
PuReTaq™ Ready-To-Go PCR Beads	see gelifesciences.com
Taq DNA Polymerase	see gelifesciences.com
illustra™ Hot Start Master Mix	25-1500-01
illustra Hot Start Mix RTG	see gelifesciences.com

Nucleotides

illustra Solution dNTPs	see gelifesciences.com
-------------------------	------------------------

mRNA and total RNA Purification

mRNA Purification Kit	27-9258-01
mRNA Purification Kit	27-9258-02
illustra RNAspin Isolation Kits	see gelifesciences.com
QuickPrep Micro mRNA Purification Kit	27-9255-01
QuickPrep mRNA Purification Kit	27-9254-01

DNA purification

ExoProStar™	see gelifesciences.com
illustra GFX™ PCR DNA and Gel Band Purification Kit	28-9034-70
MicroSpin™ S-400 HR Columns	27-5140-01
100 Base-Pair Ladder	27-4007-01

For local office contact information, visit
www.gelifesciences.com/contact

GE Healthcare UK Limited
Amersham Place
Little Chalfont, Buckinghamshire
HP7 9NA, UK

www.gelifesciences.com



imagination at work

GE, imagination at work, and GE monogram are trademarks of General Electric Company. illustra, ExoProStar, GFX, MicroSpin, PuReTaq, and Ready-To-Go are trademarks of GE Healthcare companies.

Ready to Go RT-PCR Beads: Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,789,224, 5,618,711, and 6,127,155. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

All third party trademarks are the property of their respective owners.

© 2004-2013 General Electric Company--All rights reserved. First published Dec. 2004.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Bio-Sciences AB,
Björkgatan 30,
751 84 Uppsala,
Sweden.

GE Healthcare Europe, GmbH,
Munzinger Strasse 5,
D-79111 Freiburg,
Germany.

GE Healthcare Bio-Sciences Corp.,
800 Centennial Avenue,
P.O. Box 1327,
Piscataway, NJ 08855-1327,
USA.

GE Healthcare Japan Corporation,
Sanken Bldg., 3-25-1,
Hyakunincho, Shinjuku-ku, Tokyo 169-0073,
Japan.