

## Instructions 28-9413-25 AB

# HiPrep Phenyl FF (high sub) 16/10

# HiPrep Phenyl FF (low sub) 16/10

# HiPrep Butyl FF 16/10

# HiPrep Octyl FF 16/10

## Introduction

HiPrep™ Phenyl FF (high sub) 16/10, HiPrep Phenyl FF (low sub) 16/10, HiPrep Butyl FF 16/10, and HiPrep Octyl FF 16/10 are prepacked, ready to use columns for hydrophobic interaction chromatography (HIC). The columns provide fast, preparative separations of proteins and other biomolecules based on their hydrophobic interaction with hydrophobic groups attached to the uncharged gel. See table below for column characteristics.

### Column data

Matrix	6% highly cross-linked spherical agarose (Phenyl) 4% highly cross-linked spherical agarose (Butyl, Octyl)			
Mean particle size	90 µm			
Bed volume	20 ml			
Bed height	100 mm			
i.d.	16 mm			
Column composition	Polypropylene			
Recommended flow rate*	2–10 ml/min (60–300 cm/h)			
Maximum flow rate*	10 ml/min (300 cm/h)			
Maximum pressure over the packed bed during operation, Δp†	0.15 MPa, 1.5 bar, 22 psi			
HiPrep column hardware pressure limit†	0.5 MPa, 5 bar, 73 psi			
Storage	4°C to 30°C in 20% ethanol			
	<b>Phenyl (high sub)</b>	<b>Phenyl (low sub)</b>	<b>Butyl</b>	<b>Octyl</b>
Hydrophobic ligand	Phenyl	Phenyl	Butyl	Octyl
Ligand density (µmol/ml medium)	40	20	50	5
pH stability				
short term	2–14	2–14	2–14	2–14
long term and working range	3–13	3–13	3–13	3–13

\* Water at room temperature. Flow rate is determined by  $v \bullet \leq 10$  ml/min where v=flow rate and =viscosity.

† Many chromatography systems are equipped with pressure gauges to measure the pressure at a particular point in the system, usually just after the pumps. The pressure measured here is the sum of the pre-column pressure, the pressure drop over the medium bed, and the post-column pressure. It is always higher than the pressure drop over the bed alone. We recommend keeping the pressure drop over the bed below 1.5 bar. Setting the upper limit of your pressure gauge to 1.5 bar will ensure the pump shuts down before the medium is overpressured. If necessary, post-column pressure of up to 3.5 bar can be added to the limit without exceeding the column hardware limit. To determine post-column pressure, proceed as follows:

**To avoid breaking the column, the post-column pressure must never exceed 3.5 bar.**

1. Connect a piece of tubing in place of the column.
2. Run the pump at the maximum flow you intend to use for chromatography. Use a buffer with the same viscosity as you intend to use for chromatography. Note the back pressure as total pressure.
3. Disconnect the tubing and run at the same flow rate used in step 2. Note this back pressure as pre-column pressure.
4. Calculate the post-column pressure as total pressure minus pre-column pressure. If the post-column pressure is higher than 3.5 bar, take steps to reduce it (shorten tubing, clear clogged tubing, or change flow restrictors) and perform steps 1–4 again until the post-column pressure is below 3.5 bar. When the post-column pressure is satisfactory, add the post-column pressure to 1.5 bar and set this as the upper pressure limit on the chromatography system.


## First time use

Ensure an appropriate pressure limit has been set.

Equilibrate the column for first time use or after long-term storage by running:

- a) 100 ml elution buffer (low salt) at 5 ml/min (see section "Choice of elution buffer recommendations").
- b) 100 ml of start buffer (high salt) at 5 ml/min.

These HiPrep columns can be used directly on ÄKTAdesign™ system without the need for extra connectors.



HiPrep Phenyl FF (high sub) 16/10,  
HiPrep Phenyl FF (low sub) 16/10,  
HiPrep Butyl FF 16/10 or  
HiPrep Octyl FF 16/10

### Try these conditions first

Hydrophobic interaction chromatography is usually performed with moderately high concentrations of salts in the start buffer (salt promotes adsorption) and elution is achieved by a linear or stepwise decrease in concentration of the salt.

Start buffer	0.05 M sodium phosphate buffer, 1 M (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> , pH 7.0
Elution buffer	0.05 M sodium phosphate buffer, pH 7.0
Flow rate	5 ml/min at room temperature
Gradient	0–100% elution buffer in 200 ml (10 CV)

### Equilibration before a new run

Proceed according to the instructions in section "First time use". Please read the back of these instructions for information on optimizing a separation.

## Buffers and solvent resistance

De-gas and filter all solutions through a 0.45 µm filter to increase column life-time.



### Daily use

All commonly used aqueous buffers, pH 3–13  
Guanidine hydrochloride, up to 6 M  
Urea, up to 8 M (not tested for butyl and octyl media)



### Cleaning

Sodium hydroxide, up to 1 M  
Ethanol, up to 70%  
Isopropanol up to 30%



### Avoid

Solutions <pH 2  
Phenol

## Sample preparation

Dissolve the sample in start buffer (high salt), filter through 0.45 µm or centrifuge at 10 000 × g for 10 min.



## Delivery/storage

The column is supplied in 20% ethanol. If the column is to be stored for more than two days after use, clean the column according to the procedure described under "Cleaning-in-Place (CIP)". Then equilibrate with at least 100 ml of 20% ethanol or 0.01 M NaOH at a flow rate of 5 ml/min.

**Note:** HiPrep columns cannot be opened or refilled.



## Choice of buffer

All standard aqueous buffers can be used.

When selecting salt for the start buffer refer to the Hofmeister series (see below).

Increasing the salting-out effect strengthens the hydrophobic interactions, whereas increasing the chaotropic effect weakens them.

Increasing salting-out effect										
<b>Anions:</b>	PO <sub>4</sub> <sup>3-</sup>	SO <sub>4</sub> <sup>2-</sup>	CH <sub>3</sub> COO <sup>-</sup>	Cl <sup>-</sup>	Br <sup>-</sup>	NO <sub>3</sub> <sup>-</sup>	ClO <sub>4</sub> <sup>-</sup>	I <sup>-</sup>	SCN <sup>-</sup>	
<b>Cations:</b>	NH <sub>4</sub> <sup>+</sup>	Rb <sup>+</sup>	K <sup>+</sup>	Na <sup>+</sup>	Cs <sup>+</sup>	Li <sup>+</sup>	Mg <sup>2+</sup>	Ba <sup>2+</sup>		

Increasing chaotropic effect

Table 1 lists suggested volatile buffers used in cases where the purified substance has to be freeze-dried.

Table 1. Volatile buffer systems.

pH	Substances
2.3–3.5	Pyridine/formic acid
3.0–5.0	Trimethylamine/formic acid
4.0–6.0	Trimethylamine/acetic acid
6.8–8.8	Trimethylamine/HCl
7.0–8.5	Ammonia/formic acid
8.5–10.0	Ammonia/acetic acid
7.0–12.0	Trimethylamine/CO <sub>2</sub>
8.0–9.5	Ammonium carbonate/ammonia
8.5–10.5	Ethanolamine/HCl

## Optimization

Perform your first run according to "Try these conditions first".

If the obtained results are unsatisfactory, consider the following:

Action	Effect
Change salt concentration	Higher salt concentration increases retention time Lower salt concentration decreases retention times
Increase gradient volume	May improve resolution
Decrease flow rate	Improves resolution
Increase temperature	Increases retention times
Change pH	Changes selectivity
Increase salt concentration	May increase capacity
Change to a chromatography medium with a different ligand	Selectivity change



## Cleaning-in-Place (CIP)

### Regular cleaning

Regenerate the column after each run by rinsing it with 100 ml distilled water at a flow rate of 5 ml/min at room temperature to elute material still bound to the column.

Re-equilibrate the column with at least 100 ml start buffer at a flow rate of 5 ml/min at room temperature until the UV baseline and pH/conductivity values are stable.

### More rigorous cleaning

Reverse the flow direction and run the following sequence of solutions at a flow rate of 5 ml/min at room temperature:

- 80 ml of a 1 M NaOH solution (removes precipitated proteins, hydrophobically bound proteins, and lipoproteins from the column) followed by 80 ml distilled water.
- 80 ml of 70% ethanol or 30% isopropanol (removes proteins, lipoproteins, and lipids that are strongly hydrophobically bound to the column) followed by 60 ml distilled water.

After cleaning, equilibrate the column with approximately 100 ml start buffer at a flow rate of 5 ml/min at room temperature before use.

**Note:** HiPrep columns cannot be opened or refilled.





## Troubleshooting

Symptom	Remedy
Increased back pressure	Reverse the flow direction and pump 100 ml elution buffer at a flow rate of 5 ml/min at room temperature through the column. Return to normal flow direction and run 100 ml buffer at a flow rate of 5 ml/min through the column. If back pressure is not decreased, reverse the flow direction again and follow the rigorous cleaning instructions.
Loss of resolution and/or decreased sample recovery	Follow the procedure described in the section "More rigorous cleaning".
Air in the column	Reverse the flow direction and pump 100 ml of well de-gassed start buffer through the column at a flow rate of 5 ml/min at room temperature.

## Column performance control

We recommend checking the column performance at regular intervals. Figure 1 describes how to check the function of the four different prepacked HiPrep 16/10 HIC columns.

Sample	Cytochrome C (10 mg/ml) Ribonuclease A (30 mg/ml) Lysozyme (10 mg/ml) -chymotrypsinogen (10 mg/ml)
Sample volume	2 ml
Start buffer	100 mM sodium phosphate, 1.7 M (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> , pH 7.0
Elution buffer	100 mM sodium phosphate, pH 7.0
Flow rate	2 ml/min, 60 cm/hr (room temperature)
Gradient	0–100% elution buffer in 200 ml (10 CV)
Instrumentation	ÄKTExplorer™ 100



**Fig 1.** Typical chromatogram from a function test of **a)** HiPrep Phenyl FF (high sub) 16/10, **b)** HiPrep Phenyl FF (low sub) 16/10, **c)** HiPrep Butyl FF 16/10 and **d)** HiPrep Octyl FF 16/10.

## Intended use

The HiPrep Phenyl FF (high sub) 16/10, HiPrep Phenyl FF (low sub) 16/10, HiPrep Butyl FF 16/10, and HiPrep Octyl FF 16/10 are intended for research use only, and shall not be used in any clinical or *in vitro* procedures for diagnostic purposes.

## Ordering information

Product	No. per pack	Code No.
HiPrep Phenyl FF (high sub) 16/10	1 x 20 ml	28-9365-45
HiPrep Phenyl FF (low sub) 16/10	1 x 20 ml	28-9365-46
HiPrep Butyl FF 16/10	1 x 20 ml	28-9365-47
HiPrep Octyl FF 16/10	1 x 20 ml	28-9365-48

Companion products	No. per pack	Code No.
HiTrap™ HIC Selection Kit	7 x 1 ml	28-4110-07
HiTrap Phenyl FF (high sub)	5 x 1 ml	17-1355-01
HiTrap Phenyl FF (high sub)	5 x 5 ml	17-5193-01
HiTrap Phenyl FF (low sub)	5 x 1 ml	17-1353-01
HiTrap Phenyl FF (low sub)	5 x 5 ml	17-5194-01
HiTrap Phenyl HP	5 x 1 ml	17-1351-01
HiTrap Phenyl HP	5 x 5 ml	17-5195-01
HiTrap Octyl FF	5 x 1 ml	17-1359-01
HiTrap Octyl FF	5 x 5 ml	17-5196-01
HiTrap Butyl FF	5 x 1 ml	17-1357-01
HiTrap Butyl FF	5 x 5 ml	17-5197-01
HiTrap Butyl-S FF	5 x 1 ml	17-0978-13
HiTrap Butyl-S FF	5 x 5 ml	17-0978-14
HiTrap Butyl HP	5 x 1 ml	28-4110-01
HiTrap Butyl HP	5 x 5 ml	28-4110-05
HiPrep 26/10 Desalting	1 x 53 ml	17-5087-01
HiPrep 26/10 Desalting	4 x 53 ml	17-5087-02

Accessories	No. per pack	Code No.
HiTrap/HiPrep 1/16" male connector for ÄKTAdesign	8	28-4010-81
To connect columns with 1/16" connections to FPLC™ System: Union M6 female/1/16" male	5	18-3858-01

Related literature	Code No.
Handbook, Hydrophobic Interaction Chromatography and Reversed Phase Chromatography, Principles & Methods	11-0012-69
Prepacked chromatography columns for ÄKTAdesign and Ettan™ LC systems, Selection guide	28-9317-78

## Further information

For more information, please visit  
[www.gelifesciences.com/protein-purification](http://www.gelifesciences.com/protein-purification)  
[www.gelifesciences.com/purification\\_techsupport](http://www.gelifesciences.com/purification_techsupport)

Refer also to the handbook above, contact our technical support team, or your local representative.

[www.gelifesciences.com/protein-purification](http://www.gelifesciences.com/protein-purification)  
[www.gelifesciences.com/purification\\_techsupport](http://www.gelifesciences.com/purification_techsupport)

GE Healthcare Bio-Sciences AB  
 Björkgatan 30  
 751 84 Uppsala  
 Sweden



GE Healthcare UK Limited  
 Amersham Place  
 Little Chalfont  
 Buckinghamshire, HP7 9NA  
 UK  
 GE Healthcare Europe GmbH  
 Munzinger Strasse 5  
 D-79111 Freiburg  
 Germany

GE Healthcare Bio-Sciences Corp.  
 800 Centennial Avenue, P.O. Box 1327  
 Piscataway, NJ 08855-1327  
 USA  
 GE Healthcare Bio-Sciences KK  
 Sanken Bldg., 3-25-1, Hyakunincho  
 Shinjuku-ku, Tokyo 169-0073  
 Japan

GE, imagination at work and GE monogram are trademarks of General Electric Company.  
 Ettan, HiPrep, HiTrap, ÄKTAdesign, FPLC, and Drop Design are trademarks of GE Healthcare companies.  
 All third party trademarks are the property of their respective owners.  
 © 1998–2008 General Electric Company – All rights reserved.  
 First published Feb. 1998

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.