

Amine Coupling Kit

Product description

Order code:	BR-1000-50
Contents:	<ul style="list-style-type: none">• N-Hydroxysuccinimide (NHS), 115 mg• 1-Ethyl-3-(3-dimethylaminopropyl) carbodiimide hydrochloride (EDC), 750 mg• 1 M Ethanolamine hydrochloride-NaOH pH 8.5, 10.5 mL
Storage:	2°C to 8°C
For use with:	All Biacore systems except Biacore 4000, Biacore A100 and Biacore S51.
Capacity:	The kit contains reagents sufficient for approximately 80 immobilizations.
Safety:	For use and handling of the product in a safe way, please refer to the Safety Data Sheet.

Note: *For in vitro use only.*

Note: *For Biacore 4000, Biacore A100 and Biacore S51, use Amine Coupling Kit, type 2 (BR-1006-33).*



Description

Amine Coupling Kit is designed to be used with carboxylated sensor surfaces from GE. Ligands are coupled via primary amine groups to the sensor surface for subsequent analysis.

Reagent preparation

Step	Action
1	Dissolve the EDC and NHS by adding 10.0 mL of filtered, deionized water to each vial.
2	Cap vials tightly and agitate until the solids are completely dissolved.
3	Dispense the EDC and NHS solutions separately in aliquots for storage.
4	Store the aliquots of EDC and NHS frozen at -18°C or below. Ethanolamine may be stored at 2°C to 8°C .

Note:

Use aliquots within two months.

Note: *Suitable vials for preparing aliquots for storage can be ordered from GE, for example, Plastic Vials 7 mm (BR-1002-12) and Caps (BR-1002-13), or Plastic Vials and Caps 11 mm (BR-1002-14).*

Handling

Step	Action
1	Thaw frozen aliquots and agitate them gently to make sure that the solutions are homogeneous.
2	For Biacore systems with autosampler: Place the aliquots separately in the autosampler as instructed in the software. For Biacore systems without autosampler: Mix equal volumes of EDC and NHS immediately before use.

For local office contact information, visit
www.gelifesciences.com/contact

GE Healthcare UK Limited
Amersham Place
Little Chalfont
Buckinghamshire, HP7 9NA
United Kingdom

www.gelifesciences.com/sampleprep

GE, GE monogram and Biocore are trademarks of General Electric Company.

© 2005-2016 General Electric Company.
First published Oct. 2005.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Bio-Sciences AB
Björkgatan 30, 751 84 Uppsala, Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare Bio-Sciences Corp.
100 Results Way, Marlborough, MA 01752, USA

GE Healthcare Dharmacon, Inc.
2650 Crescent Dr., Lafayette, CO 80026, USA

HyClone Laboratories, Inc.
925 W 1800 S, Logan, UT 84321, USA

GE Healthcare Japan Corporation
Sanken Bldg. 3-25-1, Hyakunincho Shinjuku-ku, Tokyo 169-0073, Japan

