

# illustra™ ExoProStar™ S

ExoProStar **S** is optimized to purify PCR and sequencing setup reactions quickly, efficiently, and reliably.

ExoProStar **S** contains illustra Shrimp Alkaline phosphatase and Exonuclease I, formulated to work together to remove unincorporated primers and nucleotides from amplification reactions in preparation for sequencing, cloning, genotyping or further DNA modification reactions.

- Enzymes optimized to work together for high efficiency removal of unincorporated primers and nucleotides
- Enzymes provided in two separate tubes, just two simple pipetting steps are needed to prepare the reaction
- Fast 15 min protocol
- Scalable for different reaction sizes
- No loss of PCR product
- Easy to automate
- Complete heat inactivation of the enzymes within 15 min




Fig 1. Schematic representation of the PCR cleanup process using ExoProStar **S**.




## Optimized for efficient primer digestion

ExoProStar **S** has been optimized for highly efficient primer digestion, helping to improve the quality of downstream analysis.

## No loss of PCR product


The use of an enzymatic digestion approach to clean up amplification reactions reduces losses of PCR product. The process has no intermediate transfer steps, spin columns or binding matrix to retain your PCR product, and double-stranded DNA is left intact by the Exonuclease I and Alkaline Phosphatase enzymes. The size of the PCR fragment does not affect the cleanup efficiency of the reaction (Fig 2).



**Fig 2.** Agarose gel electrophoresis of different size PCR products pre- and post-digestion with ExoProStar **S**. The samples were digested for 5 min at 37°C followed by denaturation at 80°C for 10 min according to kit protocol. We did not detect any loss of PCR product in any of the samples.

## High-quality sequencing results

Removal of unincorporated primers and nucleotides is essential to high-quality DNA sequencing. Failure to fully remove these components leads to high background signals and miscalling of bases. With ExoProStar **S**, Phred20 quality scores were routinely achieved at read lengths > 800 bp, equivalent to or better than other approaches to sample preparation (Fig 3).




**Fig 3.** The importance of sample clean up before DNA sequencing is illustrated in the comparison between panel A, showing PCR sequence quality following treatment with ExoProStar **S** and panel B showing sequence quality without this treatment. Read length, base calling and sequence quality are significantly improved by the use of ExoProStar **S**.

## Heat inactivation of ExoProStar **S** enzymes

Downstream operations can be adversely affected by the presence of active Exonuclease I or Alkaline Phosphatase in the PCR product following digestion. It is therefore essential that the enzymes are effectively denatured during the post-digestion heating step.

Some alkaline phosphatase enzymes are tolerant of high temperature treatment and may retain some activity causing problems in later processes. Shrimp Alkaline Phosphatase has been optimized to be quickly denatured, reducing the risk of downstream interference (Fig 4).



**Fig 4.** Temperature denaturation profile of Shrimp Alkaline Phosphatase at 65°C showing rapid and complete denaturation within 10 min. ExoProStar **S** protocol recommends denaturation of the enzyme components at 80°C, providing greater confidence in the inactivation of both enzymes prior to further downstream processes.

## Kit components and storage

ExoProStar **S** kit contains one tube of illustra Exonuclease I and one tube of Shrimp Alkaline Phosphatase. The kit is supplied on dry ice and should be stored at -20°C. Enzymes can be sub-aliquoted if required for storage convenience and should be maintained on ice during reaction setup.

## Ordering information

### ExoProStar S Enzymatic PCR and Sequence Reaction Cleanup Kit

Quantity	Code number
20 reactions	US79002
100 reactions	US79010
500 reactions	US79050
2000 reactions	US79200
5000 reactions	US79500

## Related products

### Amplification

Product	Quantity	Code number
dNTP set (100 mM each A,C,G,T)	4 × 100 µmol	28-4065-52
Illustra Ready-To-Go™ RT-PCR Beads (0.2 ml hinged tube with cap)	96 reactions	27-9259-01
Illustra PuReTaq™ Ready-To-Go PCR Beads (0.2 ml hinged tube with cap)	96 reactions	27-9559-01
Illustra Hot Start Mix RTG™ (0.2 ml tubes, 12 × 8 strip wells)	96 reactions	28-9006-53
Taq DNA Polymerase (cloned)	4 × 250 units	27-0798-05

### DNA labeling

Cy™5 dUTP	250 nmol	PA55032
Cy3 dUTP	250 nmol	PA53032
Cy5 dCTP	250 nmol	PA55031
Cy3 dCTP	250 nmol	PA53031
CyDye™ Post-Labeling Reactive Dye Pack	12 × Cy3 12 × Cy5	RPN5661

### DNA purification

illustra blood genomicPrep Mini Spin Kit	50	28-9042-64
illustra tissue and cells genomicPrep Mini Spin Kit	50	28-9042-75
illustra bacteria genomicPrep Mini Spin Kit	50	28-9042-58

### DNA cleanup

illustra GFX™ PCR DNA and Gel Band Purification Kit	100 purifications	28-9034-70
illustra GFX 96 PCR Purification Kit	10 × 96 well plates	28-9034-45
illustra MicroSpin™ S-400 HR columns	50	27-5140-01
illustra MicroSpin S-300 HR columns	50	27-5130-01

### Enzymes

illustra Shrimp Alkaline Phosphatase	500 units	E70092Y
illustra Exonuclease I	2500 units	E70073Z

For local office contact information, visit  
[www.gelifesciences.com/contact](http://www.gelifesciences.com/contact)

[www.gelifesciences.com/ExoProStarS](http://www.gelifesciences.com/ExoProStarS)

GE Healthcare UK Limited  
Amersham Place  
Little Chalfont  
Buckinghamshire, HP7 9NA  
UK



imagination at work

GE, imagination at work, and GE monogram are trademarks of General Electric Company.

Cy, CyDye, ExoProStar, GFX, illustra, ImageQuant, MicroSpin, PuReTaq, Ready-To-Go, RTG, and Typhoon are trademarks of GE Healthcare companies.

Exonuclease I and Alkaline Phosphatase method of use is covered by US patent number 5723295 in the name of GE Healthcare bio-Sciences Corp.

© 2013 General Electric Company—All rights reserved.  
First published Jul. 2013.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Bio-Sciences AB  
Björkgatan 30  
751 84 Uppsala  
Sweden

GE Healthcare Europe, GmbH  
Munzinger Strasse 5  
D-79111 Freiburg  
Germany

GE Healthcare Bio-Sciences Corp.  
800 Centennial Avenue, P.O. Box 1327  
Piscataway, NJ 08855-1327  
USA

GE Healthcare Japan Corporation  
Sanken Bldg., 3-25-1, Hyakunincho  
Shinjuku-ku, Tokyo 169-0073  
Japan